# Dynamic Mixing Chamber Dynamische Mischkammer

Manual / Handbuch

V7119-1, 03/2007

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert Knauer GmbH  
Hegauer Weg 38  
14163 Berlin, Germany  
Tel.: +49 (0)30 809 727 0  
Fax.: +49 (0)30 801 50 10  
E-Mail: [info@knauer.net](mailto:info@knauer.net)  
Internet: [www.knauer.net](http://www.knauer.net)

## CONTENTS

Using this Manual .....	3
Conventions in this manual .....	3
SOP's in this manual .....	3
General Description .....	4
Preparing the Dynamic Mixing Chamber for Operation .....	4
Unpacking .....	4
Standard Delivery List .....	4
Spare parts package for the mixing chamber .....	4
Operating the Dynamic Mixing Chamber .....	5
Description and Principle of Function .....	5
Installation .....	6
Connections: .....	6
Size of the mixing chamber .....	6
Configuration Examples of the Dynamic Mixing Chamber .....	8
Specifications of the Mixing Chamber .....	9
Warranty statement .....	18
Declaration of conformity .....	19

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## Using this Manual

This manual refers to the Dynamic Mixing Chamber.

### Conventions in this manual

Important hints are marked by the marginal hand symbol.

Special warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates helpful advice's.

### SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks for operating the Dynamic Mixing Chamber. They include step-by-step instructions assisting the user through all routine tasks during operation. They can be used for documentation purposes: They can be copied, applied, signed, and filed to document the performance of the instrument.

It is very important to follow all instructions and SOP's in this manual in order to operate the instrument and accessories. This ensures proper results and longevity of your equipment.

SOP 1	Mounting the capillaries to the mixing chamber	6
SOP 2	Pressure check of the mixing chamber	6
SOP 3	Changing the mixing chamber volume	7

## General Description

The KNAUER Dynamic Mixing Chamber was developed in order to achieve the most homogeneous mixture of different solvents with the smallest possible volumes during gradient elution.

It can be delivered in an analytical version, suitable for flow rates up to 10 ml/min and in preparative version for flow rates up to 100 ml/min. In the following the descriptions refer generally to the analytical version. Value deviations of the preparative one will be added in brackets.

**Usually the mixing chamber is made of stainless steel. For biocompatible work it also can be delivered with wetted parts made of titanium or peek.**

## Preparing the Dynamic Mixing Chamber for Operation

### Unpacking

After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transportation. If necessary, contact the carrier to claim any damages.

Use the "Standard delivery list" to check if the electric valve drive is complete. Please contact our service department if anything is missing or if you need support. It is important to fill out the guarantee card and return it to us as soon as possible.

### Standard Delivery List

1 Dynamic Mixing Chamber

Standard accessories consists of:

- 3 Stainless steel screws (M 4 x 12 mm)
- 3 Stainless steel screws (M 4 x 25 mm)
- 3 Replacement seals
- 2 Replacement sieves made of stainless steel
- 1 Hexagonal key, size 3
- 1 Operating instructions

### Spare parts package for the mixing chamber

The package consists of:

- 4 Replacement sieves (stainless steel or teflon)
- 6 Replacement Seals (stainless steel or titanium)

## Operating the Dynamic Mixing Chamber

### Description and Principle of Function


Fig. 1 *General View of the KNAUER Mixing Chamber*  
 1 Inlets    3 4-way capillary connector  
 2 Outlet    4 Power switch

The chamber of the instrument (Fig. 1) is composed of four (preparative version three) parts and is set on a housing together with a four-way 1/16" (1/8") capillary connector, which is used as inlet connection. The chamber and the connector are connected by a stainless steel capillary. The four-way capillary connector already functions as a premixer.

The analytical and preparative versions look very similarly. They only differ in size and the number of possible intermediate parts, only one for the preparative version.

The motor operates at a speed of 500 rpm which guarantees an intensive mix. Because of its overall technical concept, the mixing chamber is maintenance-free.

The mixing chamber consists of a lower and upper part as well as of two (one) intermediate pieces (see Fig. 2) and is composed of two sections. The lower half effects a fast and intensive mixing by means of the Teflon-coated stirrer. The upper half is tightly packed with glass wool and serves for additional mixing of the solvent. Without this static mixer, a fluctuating solvent composition can result from relatively strong pulsation or an uneven transport by the pumps, especially in a low flow rate range. This can cause a noisy baseline during the course of a gradient elution.

The sieves in the intermediate pieces protect the HPLC columns against solvent particles and dust from the pumps.

## Installation

**Connect the instrument to the power supply only after installation has been completed. Always disconnect the power supply before opening the instrument or the housing.**

### Connections:

#### SOP 1 Mounting the capillaries to the mixing chamber

1. Connect the pump(s) to the inlets of the 4-way connector with capillaries using bushings and ferrules. First, tighten the bushings with which the pump(s) are connected to the four-way connector by hand, then 1/2 -3/4 of a turn with a wrench.

**For high pressure gradients, up to three pumps can be directly connected to the four-way connector on the mixing chamber. To do this, the three ports on the connector equipped with blind fittings are to be used. For low pressure gradients where only one pump is required, any one of the three inlets can be used. The remaining two inlets must be closed with blind fittings.**

2. Connect the outlet of the mixing chamber (2, Fig. 1) to the injection valve, then proceed as above.
3. Now connect the instrument to a grounded power supply outlet.
4. Before using the mixing chamber for analyses, rinse it with approx. 100 ml methanol or i-propanol. Check the chemical compatibility of the eluents to be used with the materials in the dynamic mixing chamber.

#### SOP 2 Pressure check of the mixing chamber

1. Set the following parameters at your pump:

Flow rate:	50 (250) µl/min
Maximum pressure:	250 bar with 10 ml pump head 50 bar with 50 ml pump head
2. Connect the pump head outlet to one of the mixing chamber inlets. Blind the other inlet and the outlet of the mixing chamber.
3. Start the eluent delivery and wait for the pressure stop of the pump. The displayed pressure shall not drop more than 5% within three minutes.
4. Inspect during this time all connections for possible leakages.

**Even the smallest leakage makes a correction necessary. Repeat in this case the pressure check.**

### Size of the mixing chamber

The mixing chamber is assembled as a complete unit for delivery so that it has the largest possible volume. It is therefore suitable for use in the semi-preparative range. Should very short response times be required, remove one or both of the intermediate pieces according the following SOP 3. Removal of the static mixing section (upper intermediate piece filled with glass wool) should scarcely be necessary since, although this would further shorten the response time, it does not affect the mixing speed (see Configuration Examples of the Dynamic Mixing Chamber).

**SOP 3 Changing the mixing chamber volume**

1. Disconnect the power supply of your mixing chamber.
2. Disconnect the capillary of the chamber outlet.
3. Loosen and remove the three mounting screws.
4. Remove the chamber upper piece.
5. According to your requirements remove or set in the chamber intermediate piece(s).
6. If all intermediate pieces are removed place a sieve (7) directly on the lower part.
7. Set in the chamber upper piece and fix it with the three mounting screws. If the intermediate piece is removed, the three shorter screws out of the accessory kit must be used..
8. Connect the outlet capillary to the chamber outlet.
9. Perform the SOP 2 to ensure the mounting without any leakage.


Fig. 2 Construction of the mixing chamber

- |                                       |                           |
|---------------------------------------|---------------------------|
| 1 lower part                          | 5 Teflon seal             |
| 2 Intermediate part (dynamic section) | 6 Stirrer (magnetic core) |
| 3 Intermediate part (static section)  | 7 Sieve                   |
| 4 upper part                          | 8 Glass wool filling      |

## Configuration Examples of the Dynamic Mixing Chamber

The role of the chambers is as follows:

- The lower chamber range effects rapid, dynamic mixing by means of a Teflon coated stirring magnet.
- The upper chamber provides additional static mixing of the solvent. This static mixer helps to minimize the pulsation and/or uneven delivery of solvent from the pumps, which can occur at lower flow rates and lead to noisy baselines, especially when a gradient is employed.

Some examples shall illustrate the possibilities for use of the dynamic mixing chamber in various configurations. They demonstrate how the actual composition of the eluent at the column inlet depends on the total mixing area volume before the column, as well as on the dead volume of the static (with glass wool) or dynamic (without glass wool) mixing chamber.

The examples are based on a system of eluent A = pure water and eluent B = methanol/water 80/20 v/v. At time zero always a switch from 100% A to 100% B is performed. The detector was directly connected to the dynamic mixing chamber. Thus, the illustrations represent the composition of the eluent at the column head.

Fig. 3 symbolizes all possible configurations of the mixing chamber. The preparative version can only be configured in the variants A, B, and C.


Fig. 3 Mixing chamber configurations

Firstly, the influence of the flow rate at a given mixing chamber volume shall be shown (Fig. 4). The higher the flow rate the steeper is the slope of the eluent gradient, limited of course by the capacity (volume) of the mixing chamber.


Fig. 4 Composition of eluent at 1.0, 0.5, and 0.2 ml/min in configuration A 

Very similar is the influence of the mixing chamber volume at a given flow rate. The smaller the mixing chamber volume the steeper is the slope of the eluent gradient.


Fig. 5 Composition of eluent in configurations **A** , **B** , and **D**  at 1 ml/min

The use of filled intermediate pieces enlarges the static mixing area in addition to the dynamic mixing volume of the lower part. The influence on the response time is obviously whereas on the curve form it is much lower.


Fig. 6 Composition of eluent in configurations **A** , **C** , and **G**  at 1 ml/min

Fig. 7 illustrates that for all possible combinations with two intermediate pieces the curve shape is mainly dependent on the lower one if it is filled or not.


Fig. 7 Composition of eluent in configurations **G** , **F** , **E** , and **D** 

Flow rate and volume of the mixing chamber affect the form of the gradient. Mixing chambers from different manufacturers can therefore produce different gradients at the start of the column. The usual details given in the literature for the gradient curve of the controlling instruments are generally not sufficient for reproduction of the same results with other mixing chambers.

## Specifications of the Mixing Chamber

	Analytical	Preparative
Chamber Volume		
Total	1740 $\mu$ l	5.9 ml
Lower part	260 $\mu$ l	1.1 ml
Dynamic intermediate	880 $\mu$ l	4.8 ml
Static intermediate	500 $\mu$ l	-
Maximum Pressure	420 bar	250 bar
Power Supply	110-220 V AC	110-220 V AC
Dimensions (WxHxD)	80 x 100 x 120 mm	80 x 130 x 125 mm
Weight	1.2 kg	1.6 kg

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## Warranty statement

The guarantee period of the Dynamic Mixing Chamber is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Claims under this guarantee are valid only if the enclosed guarantee card is returned to us at the address shown below within 14 days of receipt of the instrument. Defective Dynamic Mixing Chambers should be sent to the manufacturer for repair.

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert KNAUER GmbH  
Hegauer Weg 38  
D-14163 Berlin  
Tel: 030 – 809 727 – 0  
Fax: 030 – 801 50 10  
e-mail: [info@KNAUER.net](mailto:info@KNAUER.net)  
[www.KNAUER.net](http://www.KNAUER.net)

If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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## Declaration of conformity

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**Manufacturer's name and address:****Herstellername und -adresse**

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Order numbers, Bestellnummern

**A 0275, A0285, A0581 and A1174**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)  
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)  
EN 6100-3-2 (10/2006)  
EN 61326-1 (10/2006)

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
The product was tested in a typical configuration.

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Berlin, 2007-04-02

Alexander Bünz (Managing Director)

The CE Shield is attached to the rear of the instrument.

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